



MULTIPLE ORTHOTIC OPTIONS

The New Contour™ LSO by Aspen provides a simple, effective, comfortable and economic spinal bracing system. Created by combining the QuikDraw™ RAP and the Contour™ Back Panel, the Contour™ LSO provides spinal bracing for patients of any size, shape or medical condition.

Simple

The Contour™ LSO is Simple to apply and use. Only one measurement is needed to size the patient. Once a size is chosen, simply overlap the ends of the QuikDraw™ and then compress by pulling the tabs using the QuikDraw™ SlickTrack™ tightening system.

Effective

The Contour™ LSO is Effective, providing a rigid layer of support for the patient. Depending on the motion restriction desired, choose from 3 back panel sizes: 18", 15" or 12" in height.

Comfortable

The Contour™ LSO is Comfortable. By providing an effective brace that is also comfortable, patient compliance is increased. This, in turn, will enhance the progression of care.

Economic

The Contour™ LSO is Economic, providing multiple braces in one spinal bracing system. This minimizes stocking requirements, keeping costs low.

L0637 Code Approved


 **Aspen**
MEDICAL PRODUCTS INC. 6481 Oak Canyon
Irvine, CA 92618
800-295-2776
www.aspenmp.com

Innovative Patient Care



PRODUCT INFORMATION


CONTOUR™ BACK PANEL



CATALOG NUMBER

	WHITE	BLACK
18" Back Panel Height = 18"	992118	992218
15" Back Panel Height = 15"	992115	992215
12" Back Panel Height = 12"	992112	992212
18" Replacement Pad Set	992128	
15" Replacement Pad Set	992125	
12" Replacement Pad Set	992122	


QUIKDRAW™ PRO




CATALOG NUMBER

SIZE	WHITE	BLACK
QUIKDRAW™ EXTRA SMALL	993080	993081
QUIKDRAW™ SMALL	993010	993011
QUIKDRAW™ MEDIUM	993020	993021
QUIKDRAW™ LARGE	993030	993031
QUIKDRAW™ EXTRA LARGE	993040	993041
QUIKDRAW™ EXTRA EXTRA LARGE	993090	993091
COLD THERAPY GEL PACK	993005	
EXTENSION PANEL	993060	993061


Rigid Anterior Panel (RAP)



Regular RAP White
Cat. No. 993100
Regular RAP Black
Cat. No. 993200
Regular RAP Pad
Cat. No. 993110



Large RAP White
Cat. No. 993101
Large RAP Black
Cat. No. 993201
Regular RAP Pad
Cat. No. 993110



Extra Large RAP
Cat. No. 993102
Extra Large RAP Pad
Cat. No. 993112

QuikDraw™ + RAP (R, L, XL) + Contour™ Back Panel = Contour™ LSO

	EXTRA SMALL	SMALL	MEDIUM	LARGE	EXTRA LARGE	EXTRA EXTRA LARGE
MEN'S PANT SIZE	20 to 24 in.	24 to 28 in.	28 to 32 in.	32 to 36 in.	36 to 42 in.	42 to 54 in.
WOMEN'S PANT SIZE	Child to 0	1 to 3	5 to 9	11 to 15	17 to 21	23 to 33